

Blood Bank Refrigerators

B Medical Systems | B Range

For the safe storage of blood and blood components being subject to cold chain and temperature sensitivity > according to DIN 58371 | ÖNORM K 2030













SAVING LIVES
THROUGH RELIABLE
AND INNOVATIVE
TECHNOLOGY

Blood Management Solutions



Safety Standards **B Medical Systems**











The Safety Standards developed by B Medical Systems define certain significant technical features of a product. These ensure the safe storage of the preparations as well as setting the highest standards of safety for the user.

PRECISION LINE	PREMIUM LINE	
		B Medical Systems Electronics
		B Medical Systems digital display
		B Medical Systems 4,3" digital color display with capacitive buttons
	0	B Medical Systems 7" full touchscreen display
		Safety door lock (with 2 keys) and key-operated power switch ON/OFF (with 2 keys)
		Power indicator light and digital temperature indicator (display: 0.1 digits)
		Controlled fan cooling system for constant temperature and even temperature distribution across the entire refrigerating chamber
		Automatic switch-off of the evaporator fans when the door opens
		Self-contained alarm system with integrated battery takes over the alarm function and temperature value measurements in case of system failure for at least 48 hours
		Acoustic and visual alarm signal in case of temperature alarm and system failure
		The alarm history function on the electronic stores all the relevant values during a temperature alarm, such as: min., max. and average temperature [+ alarm duration for Premium Line models]
		Remote transmission alarm signal (via potential-free contact) in case of temperature alarm (change-over contact)
		Safety thermostat prevents dropping of the temperature below +2°C
		Door opening alarm (visual / acoustic)
		Designed and tested for climatic class SN (ambient temperature range +10°C to +32°C)
		Designed and tested for climatic class SN / T (ambient temperature range +10°C to +43°C)
		Interior made from stainless steel
		Internal LED lighting
		Additional remote transmission alarm signal (via potential-free contact) in case of system failure (change-over contact)
		Smooth castors for optimum flexibility of movement
	0	Ambient temperature sensor
0		B Connected - Universal software for the monitoring of refrigeration devices, including the acquisition, recording and visualization of temperature data
		DCU (Data Communication Unit) - Hardware module monitoring all operating conditions and passing them through to a central database – via local network (in combination with °B Connected)
		Ethernet interface for the visualization of all operating and control functions (hardware and software settings) via B Connected monitoring software on a peripheral device (computer)

MODELS B 51-131-291-381 & BF 261 | PRECISION LINE

Blood Bank Refrigerators

B Medical Systems | B Range

9 models • Volume 52 > 895 L • Set temperature +4°C • Climate class SN | SN/T • According to DIN 58371 | ÖNORM K 2030

In conformity with national and international guidelines, regulations for Medical Devices offering reliability, efficiency and safety at an optimal price.

Integrated multi-function monitoring electronics

Integrated B Medical Systems multi-function electronics; with easy access in the door level; offers alarms, central alarm system, set point security, battery backup, full display and many other control features; compatible with °B Connected monitored solution for data monitoring and recording.



• Lifetime warranty corrosion-free cabinet

- One piece material cabinet with a lifetime warranty for zero corrosion, highest quality polypropylene for better insulation and no sharp edges to provide easier maintenance and cleaning
- · Smart ergonomy Easier and safer handling because of heavier components placed at the bottom



High storage capacity

High storage per square meter in terms of net volume, with a very versatile inner volume utility due to multiple user friendly rails allowing modularity.




Exclusive integrated electronics handle bar (Patent pending)

Blood bank refrigerators for the safe storage of blood

and blood components at temperatures ranging from +2°C to +6°C with multiple controls, warning and data recording safety features developed as a result of over 35 years of expertise in refrigeration technology.

All functionalities are easily accessible:







High quality materials

New high quality coating, certified medical devices quality and antibacterial, high quality steel for better longevity and easy hygiene



Versatile modularity

Excellent storage capacity and modularity - Large choice of shelves and drawers offering modularity for every need.



4-layer (Precision Line) and 5-layer (Premium Line) PMMA organic glass offering highest insulation, combined with reduced weight reducing efforts to close door and constraints on hinges, offering better longevity.







MODELS B 401-501-701-901 | PREMIUM LINE





Technical DataGeneral features











	В	51	В	131	В	291	В	381	BF	261	Refrigerator 8 combined	+ Freezer
Gross / Net volume (I)	52 / 45		121 / 106	121 / 106		297 / 269		362 / 329		06	F 121 / 106	
torage capacity (blood bags)	24 (450ml)		51 (450ml)		159 (450ml)		213 (450ml)		51 (450ml)		64 (450ml)	
	31 (350ml)		79 (350ml)		235 (350ml)		313 (350ml)		79 (350ml)		104 (350ml)	
et temperature (preset)					+	4°C					-32°C	
emperature cold / warm alarm limit					+2°C	:/+6°C					-37°C / -27°	°C
fold over time (+4°C to +10°C)	2 h 30		2 h 42		2 h 12		1 h 52	1 h 52			42 min.	
Climate class (ambient temperature range)					SN (+10°	C to +32°C)						
Defrosting technique	Natural									Manual (for Freezer)		
Refrigerant type					R	500a					R290	
xternal dimensions H x W x D (mm)	670 x 495 x	575	830 x 595 x 695		1670 x 595 x 695		1980 x 595 x 695		1660 x 595 x 695			
nner dimensions H x W x D (mm)	470 x 375 x	350	630 x 475 x 470		1470 x 475 x 470		1790 x 475 x 470		630 x 475 x 470		630 x 475 x 470	
let weight with standard equipment (kg)												
Supply voltage (V)	220-240	115-127	220-240	115-127	220-240	115-127	220-240	115-127	220-240	115-127	220-240	115-127
Frequency (Hz)	50 / 60	60	50 / 60	60	50 / 60	60	50 / 60	60	50 / 60	60	50 / 60	60
Power (W)	42	46	110	82	120	85	135	94	110	82	168	209
Energy consumption (kWh/24h)	0.42 / 0.44	0.50	0.85 / 0.86	0.68	1.13 / 1.21	0.99	1.22 / 1.28	1.06	0.85 / 0.86	0.68	1.35 / 1.28	1.98
Heat emission (Kcal/h)	12	12	27	24	34	34	41	34	27	24	53	65
Compressor running time (%)	32	31	28	34	33	46	35	42	28	34	37	36
Noise level (dB(A)) (at 1m height & 1m distance)	35 / 36	37	35 / 36	38	38 / 40	39	40	38	35 / 36	38	45 / 47	47



Technical DataGeneral features









		В	401	В	501	В	701	В	901		
Gross / Net volume (I)		451 / 384		598 / 510		747 / 636		895 / 763			
Storage capacity (blood bags)		255 (450ml) 345 (350ml)		340 (450ml) 460 (350ml)		425 (450ml) 575 (350ml)					
Set	temperature (preset)	+4°C									
Ter	nperature cold / warm alarm limit				+2°	°C / +6°C					
Но	ld over time (+4°C to +10°C)										
Clii	mate class (ambient temperature range)				SN / T (+	10°C to +43°C	:)				
De	frosting technique	Natural									
Refrigerant type						R600a					
Ext	ernal dimensions H x W x D (mm)	1988 x 697 x 1035		1988 x 850 x 1035		1988 x 982 x 1035		1988 x 1140 x 1035			
Inn	er dimensions H x W x D (mm)										
Ne	t weight with standard equipment (kg)										
	Supply voltage (V)	220-240	115-127	220-240	115-127	220-240	115-127	220-240	115-12		
ver	Frequency (Hz)	50 / 60	60	50 / 60	60	50 / 60	60	50 / 60	60		
y Sa	Power (W)										
Energy Save	Energy consumption (kWh/24h)										
	Heat emission (Kcal/h)										
	Compressor running time (%)										
	Noise level (dB(A)) (at 1m height & 1m distance)										



Technical Data

Specifications

		B 51-131 & BF 261	B 291-381	B 401-501-701-901
Regulation se Precision (from -8	ensor 80°C to +180°C) ± 0.3°C	NTC		
Display senso Precision (from -8 in reference body corning 200-5CST	80°C to +180°C) ± 0.3°C v / 100 ml DOW	NTC		
Accu data / function time of the control panel when system failure		12 V - 2.3 AH	/ 48h	
Relative hum	nidity at +32°C		≤	75%
Door insulation	on (polyurethane)	80 mm PU with 4-Layer Pl	MMA organic glass	
Casing insula	ntion (polyurethane)	60 mm P	U	
	Inner / Outer body	Polypropylene (cop	oolymer UV)	
	Inner door	Polypropylene (copolymer UV)	Aluminium	
Material	Outer door	Polypropylene (cop	oolymer UV)	
	Drawer	Polycarbonate, tr	ansparent	Stainless steel (V2A - 1.4301)
	Shelf	-		Wire DIN177, PA11 coated
European Me	edical Device Directive	93 / 42 / EEC, Class IIa (m	odels: 220-240 V - 50/60 Hz)	
EMC directive	е	2014 / 30 / CE (models:	220-240 V - 50/60 Hz)	
Low voltage	directive	2014 / 35 / CE (models:	220-240 V - 50/60 Hz)	
FDA regulation	on	Class II medical devices acc. to 21CFR I	Part 864.9700 (models: 115-127 V - 60 Hz)	
cCSAus certifi	icate	Yes (models: 115-12	7 V - 60 Hz)	
GMP - Clean	room classification	Class ISO 6 / EC	C GMP B	




Equipment Standard & optional



DRAWER 3 versions: - top - with Front Cover - bottom





WIRE SHELF




		B 51	B 131	B 291	B 381	BF 261	B 401	B 501	B 701	B 90°
B Medical	Digital display	•	•	•	•	2	-	-	-	-
Systems Electronics	4.3" digital color display	-	-	-	-	-	•	•	•	•
	7" full touchscreen display	-	-	-	-	-	\circ	0	0	0
	Top Bottom Drawer	- • 1	- • 1	- • 1	- • 1	●1 ●2	-	-	-	-
	Drawer with Front Cover	1	1	5	6	2	-	-	-	-
nterior quipment	ST-Drawer with Front Cover	-	-	-	-	-	3	3	3	3
quipinent	Wire Shelf	-	-	-	-	-	2	2	2	2
	ST-Shelf	-	-	-	-	-	\circ	0	0	0
Internal LED lighting		•	•	•	•	•	•	•	•	•
B Connected -	Monitoring Software	0	0	0	0	0	•	•	•	•
CU - Data Con combination wit	mmunication Unit th °B Connected)	0	0	0	0	0	-	-	-	-
thernet interf	ace	-	-	-	-	-	•	•		
mbient temp	erature sensor	-	-	-	-	-	\bigcirc	\circ	\circ	\circ
Potential-free (a case of system fa	alarm contact _{illure}	•	•	•	•	•	•	•	•	•
Integrated port for external sensor (installed by customer)		•	•	•	•	2	•	•	•	•
Rollers Smooth castors with stabilizers		- -	2 -	2 -	2 -	2 -	- • 4	- • 4	- • 4	- • 4
Door hinge Right Left										
Wooden packaging for ocean transport / export		0	0	0	0	0	0	0	0	0







S-RACK



Accessories Standard & optional

	B 51	B 131	B 291	B 381	BF 261	B 401	B 501	B 701	B 901
Dividers Set (H 100 mm) to facilitate the handling and storage management	• 1	1	6 5	6	2	-	-	-	-
S-Rack (H 100 x W 113 x D 638 mm)	-	-	-	-	-	0	0	0	0
Temperature chart recorder pre mounted (in the form of a circular chart recorder)	0	0	0	0	0	0	0	0	0

Full functionalities at a glance

B Medical Systems | Electronics

PRECISION LINE | Multifunctional electronics with digital display and easy access, in the door handle; compatible with °B Connected monitored solution.



• Key power switch • O = OFF / I = ON

- To switch the unit on, turn the key switch to position "I". The green power LED lights up
- After the self-check, the temperature inside the device is displayed




Digital display

During normal operation the display shows:

- Left-hand side: current time and internal temperature displayed in °C (Celsius) or °F (Fahrenheit) in increments of 0.1°
 - Right-hand side: the maximum and minimum inside temperature



Function buttons





 ENTER button: use the ENTER-key to confirm the actual selection. This applies the parameter value. The change of a parameter only becomes effective when confirmed



 MUTE / BACK button: in case of an alarm, the acoustic alarm can be muted for a specified time. In Parameter-mode the same key acts a "BACK-button" to go one step back



 LIGHT button: allows to switch on or off the interior light of the device without opening the door





When an alarm occurs, the red alarm LED lights up and an acoustic signal sounds. The corresponding alarm message appears on the display alternately with the measured inside temperature. If several alarms occur simultaneously, the messages will be shown alternately.



A USB port is available for data export only. Measured values of the display sensor are stored for at least 30 days (ring buffer) together with a time stamp.



3.8°C

HAL

The new B Medical Systems Electronics

offers a wide range of adjustment and diagnostic facilities as well as additional protection / warning operations (via external alarm operations, histories and individual display signals).

PREMIUM LINE | Multi-function electronics with 4.3" digital color display and capacitive buttons, integrated at optimal level in the door handle with full commands and connections.



OFF / ON

- To switch on press the ON/OFF button for 1 second. The green LED will light up
 - Switch off protected by password.
 The door will automatically unlock before powering down



O Digital color 4,3" display

During normal operation the display shows:

- At the top: current voltage and current time
 - Left-hand side: internal temperature displayed in °C (Celsius) or °F (Fahrenheit) in increments of 0.1°
- Right-hand side: temperature of the regulation sensor. Displaying this temperature can be disabled through the settings



Acoustic and visual alarm signal

When an alarm occurs, the red alarm LED lights up and an acoustic signal sounds. The corresponding alarm message appears on the display together with the measured inside temperature. If several alarms occur simultaneously, the messages will be shown alternately.



• Capacitive buttons

- ARROW buttons: select screen items with the arrow buttons. In parameter editing mode, the arrow buttons are used to change values
- ENTER button: activates the currently selected item
- BACK button: goes one level upwards in the menu tree structure. It also serves as "escape" button
- DOOR UNLOCK button: unlocks the door of the device
- MENU button: opens the menu screen
- ALARM MUTE button: deactivates the acoustic alarm for a predefined period of time
- LIGHT button: allows to switch on or off the interior light of the device without opening the door

















Discover the ultimate control experience with the new B Medical Systems 7" full touchscreen display



Available as option

More confortable, more functionalities, no compromise!

Optimum control and protection

B Medical Systems | °B Connected & DCU

°B CONNECTED | MONITORING SOFTWARE

Universal software for the monitoring of refrigeration devices, including the acquisition, recording and visualization of temperature data.



- Unique monitoring software for the full range of Blood Management Solutions and Medical Refrigeration products
- Web-based interfaces for computers and mobile devices
- Modern design for simple and intuitive use
- Graphic display of temperature curves
- Integrated event and activity history of appliances' components
- Data recording on centralized database for long-term archiving
- Easy setting of specific alarm forwarding, via email or SMS
- Generation of reports compiling crucial data and activities
- Temperature and detailed device data export for third-party software
- Important price advantage compared to a traditional circular chart recorder and its spare parts
- REST API to access raw data directly from the database in read only mode

KEY BENEFITS:

- Unique interface for the monitoring of the full range of refrigeration devices
- Centralized database providing data access to entire customer network









- Optimised design for mobile and touch compatible devices Support of multiple
- Real time data: Current temperatures, set point, upper and lower limit, door state, compressor runtime
- · Log data: Events, alarms within date/time, event type, description of the event. Event and alarm confirmation function
- Batch management: Read in of plasma bags with a barcode scanner
 - Movement state and door control of Agitator devices
- Parameter settings of each device: Factory settings, real time clock, language, firmware version, RS485 address, IP address

> Comprehensive applications and diagnostic possibilities







DCU | DATA COMMUNICATION UNIT (PRECISION LINE)

Hardware module monitoring all operating conditions and passing them through to a central database – via local network.



- Interface connection of B Medical Systems appliances to an existing network
- Graphical user interface displaying temperature and configurating devices and alarms easily
- Direct connection to Ethernet and serial bus RS485
- Digital in/out (programmable and customer-specific use)
- Recording and storing appliance relevant data
- Integrated USB port allowing to export stored data
- The DCU combined with the °B Connected software replaces the paper temperature recorder
- All data are recorded and save in the internal storage of the DCU and backed up in °B Connected data base if connected
- Several additional self-sufficient temperature sensors (up to 4 PT1000) may be connected
- Humidity sensor input (4-20 mA)

KEY BENEFITS:

- Central system collecting relevant temperature data of the appliances and their respective operating conditions
- A number of connection facilities allowing flexible upgrades for individual projects







State-of-the-art technology for the exacting needs of the medical world

Reliable and innovative solutions for safe vaccination around the world







After Sales support and service

We strive to provide you with the highest standards of service; not only through our selected distributors and partners for all your maintenance and service but also our second line trouble shooting and after sales service. This Factory based group of engineers is there to help our partners and yourself to get the best solution for your cold storage needs.



Our Global Expertise



SAVING LIVES THROUGH RELIABLE AND INNOVATIVE TECHNOLOGY

B Medical Systems (formerly Dometic Medical Systems) has more than 35 years' experience in the medical refrigeration sector.

The company, formerly known as Electrolux Medical Systems, was founded in 1979 when the World Health Organization approached Electrolux in Vianden, Luxembourg, to create a solution for the safe storage and transport of vaccines around the world. In 2001, Electrolux Medical Systems became part of the Dometic Group, and was renamed Dometic Medical Systems. Having established a legitimate reputation in the medical equipment industry, the company has also become a global leader in vaccine cold chain.

B Medical Systems S.à r.l.

17, op der Hei L - 9809 Hosingen, Luxembourg

Tel.: (+352) 92 07 31-1 Fax: (+352) 92 07 31-300 info@bmedicalsystems.com



Luxembourg, in the heart of Europe



www.bmedicalsystems.com