





The perfect lounge ...





Easy to carry on rolls



■ **Portable:** The donor lounge is foldable.

■ Shock position: A collapsed donor can be set

into a shock position supported

by a gas tube.

■ **Donor slip:** Easy to clean and desinfect.

The material has a net structure and is comfortable

for the donor.

■ **Head rest:** The upholstery of the lounge is

integrated in the unit.

■ **Arm rest:** The armrest is adjustable on

both sides.





Shock positionGas pressure supported



Lmb Technologie GmbH

Möslstr. 17

D-85445 Schwaig

Tel: +49-8122-88096-0 Fax:+49-8122-88096-60

www.lmb.de info@lmb.de

Technical data	
Material:	Stainless steel (round pipe) to prevent scratches and accidents
Weight without arm rest:	17,8 kg
Arm rest:	1,7 kg
Max. net load:	130 kg
Dimensions:	125 x 55 x 22 cm (HxWxD)
High of seat:	58 cm
High of arm rest:	75 cm
Seat extent:	55 cm
Angle of seating:	approximately 140°